

Specification Sheet

Product Name : Tactimab

Cat. No. : 3121, 3122

General name: Affinity purified Anti Human CD3 Mouse Monoclonal Antibody

Isotype: Mouse IgG2a, Kappa

Purification: Protein G affinity chromatography

Size : 1 ml/vial containing 1 or 10 mg/ml of Tactimab

Source: Animal protein-free culture of Mouse OKT3 clone 1B6

Storage Condition: Below -20°C

Shipping Condition: -20°C

EXAMINATION		SPECIFICATION
APPEARANCE (Turbidity)		Clear
APPEARANCE (Form)		Color less clear solution* ²
pH		6.5~7.5
Purity Test		>98% by SDS Polyacrylamide Gel Electrophoresis
Protein		Determine by Absorbance at UV 280nm
Endotoxin		<0.1EU/mg by LAL test
Sterility Test	Bacteria	Negative
	Fungi	Negative
Cell culture test * ¹		PASS

*1 Human T Lymphocyte proliferation assay.

*2 Phosphate buffer solution containing 0.02% Tween 80.

Note: Research purposes only. Not for diagnostic or In vivo use.