

Specification Sheet

Product: EALI515-NK (Code# 1602, 1622, 1631, 1632)

Examination	Methods	Specification
Color	Visual test	Orange – Red
Turbidity	Visual test	Clear
pH (After equilibration to 5% CO ₂ at 37°C)	JP* standard	7.4–7.7
Osmotic pressure	Freezing point depression method	290±10 mOsm/kg H ₂ O
Endotoxin	Kinetic–chromogenic technique by JP* standard	< 0.1EU/ml
Sterility	Sterility test by JP* standard	Negative
Mycoplasma	Nucleic acid amplification test	Negative
Cell culture	Cell proliferation test by HPBL**	Company Standard

* Japanese Pharmacopoeia

** Human peripheral blood lymphocyte